

## Human Chromosome Paint Probe

### Needed But Not Provided

- 
- Tween 20
  - Distilled water
  - Formamide
  - 20x SSC
  - 80%, 100% ethanol (stored at room temperature)  
Prepare by adding dH<sub>2</sub>O to 800ml pure ethanol to a final volume of 1L.
  - 70% Ethanol stored in -20°C
- 

### Reagents preparation:

#### Day 1

##### Denaturation solution (70% formamide /2SSC)

Add 35ml formamide, 10ml distilled H<sub>2</sub>O, 5ml 20xSSC Adjust pH to 7.0 using HCL, heat to 72°C.

#### Day 2

##### Rapid wash (0.4 X SSC solutions)

Add 1ml 20X SSC, 49ml distilled water. Mix well and heat to 74°C.

##### Washing solution II (4 X SSC/0.1%Tween 20)

Add 100ml 20X SSC, 400ml distilled water, 0.5ml Tween 20.

## Protocol

#### Day 1

##### A) Chromosome denaturation

1. Put the slides in 2XSSC at RT for 2 min and then dehydrate in Ethanol series: 70%, 80% and 100%, 2 min. each. Air dry.
2. Heat 40ml of denaturation solution to 70°C ( $\pm 2^\circ\text{C}$ ) in a glass Coplin jar. Place slides in the solution for 1.5 minutes. DO NOT OVER DENATURE, some samples denature in 60 seconds. Hot plate can also be used for denaturation: put 100ul of the denaturation solution on the slide, cover with a cover glass and put on a slide warmer at 72°C ( $\pm 2^\circ\text{C}$ ) for 1.5 minutes.
3. Immediately place slides in Cold 70%, and in 80% and 100% ethanol, 2 min. each. Air dry.

##### B) Probe denaturation and hybridization

1. Centrifuge briefly the content of the probe mixture, take 10ul for each slide and denature the probe by incubation at 80°C in a water bath for 7 minutes.
2. Put in a water bath at 37°C for 10 minutes.
3. Add 10ul from the denature probe mixture to the denatured chromosome preparation.

***For Research Use Only***

4. Place an 18 x 18mm<sup>2</sup> cover slip over the probe mix, being careful not to trap air bubbles under the cover slip. Seal the edges with rubber cement. Transfer the slide to a humidified chamber or container and place in incubator or baking oven set at 37 °C for 12-16 hours.

**Alternatively: Co-denaturation can be used: apply 10 ul from the probe, put a cover glass (18X18mm) and seal with rubber cement. Denature sample and probe together on a hot plate at 74 °C for 4 minutes. Place in an incubator or baking oven set at 37 °C for 12-16 hours.**

## Day 2

### Detection

**Note: During the whole procedure the slides should remain wet and protected from direct light.**

1. Remove slides from the humidified chamber and carefully remove the rubber cement.
2. Transfer the slides to a Coplin jar containing 0.4XSSC. Wash slides in 0.4XSSC at 74 °C ( $\pm$  2°C) for 3-5 min. Dip slides in washing solution II (4XSSC/ 0.1% Tween 20) for 2 minutes.
3. Put 20ul of antifade solution with DAPI place a cover glass (24X60mm<sup>2</sup>) over the surface. Try to remove any air bubbles that may have formed.

***For Research Use Only***